JEILMEDICAL

CE  
1639

# Dual Top<sup>®</sup> ANCHOR SYSTEM

## Instruction For Use

- Dual Top Anchor System
- Sterile Dual Top Anchor System

---

Issue Date: 2023. 04

Print Date:

#497-41001-08

## Table of Content

English .....	03
Deutsch .....	10
Italiano .....	18
FRANÇAIS .....	25
Español .....	32
Čeština .....	39
Ελληνικά .....	46
Magyar .....	53
Nederlands .....	60
Polski .....	68
Português .....	75
Română .....	82
Svenska .....	89
Türkçe .....	96
Slovensko .....	103

**[DESCRIPTION]**

The Dual Top Anchor System was designed to provide a firm anchorage in the oral cavity for an orthodontic treatment. Plates are manufactured from titanium (ASTM F67) and Screw are manufactured from titanium alloy (ASTM F136). The Dual Top Anchor System is designed to apply to various orthodontic tools and diverse enough to satisfy various clinical cases.

**[DIRECTIONS FOR USE]****1. Dual Top Screw****1-1. Screw insertion without drilling**

- ① Select a screw and a driver.
- ② Administer a topical anesthetic to a graft site.
- ③ Attach the selected screw to the driver.
- ④ Place the screw straight on the mucosa, and then turn the handle as pressing it so that the tip of the screw can dig a hole in the surface cortical bone. (While inserting the screw, turn it slowly to prevent the driver axis from shaking.)
- ⑤ Stop inserting when the head at the lower end of the screw meets the mucosa, and remove the driver from the screw.

**1-2. Screw insertion after drilling**

- ①, ② The same as above.
- ③ Attach pilot drill to a hand piece, and start drilling. (1,300 rpm or less. The drilling area should be cooled by physiologic saline.)
- ④ Attach the selected screw to the driver.
- ⑤ Turn the handle, fitting the hole in the cortical bone. (While inserting the screw, turn it slowly to prevent the driver axis from shaking.)
- ⑥ Stop inserting when the head at the lower end of the screw meets the mucosa, and remove the driver from the screw.

**2. Anchor Plate**

- ① After the diagnosis was concluded about the orthodontic treatment, the surgeon selects the surgical site.
- ② The device set on the selected site and then fixate the plate by using the screw(20-AT Series). In some case the surgeon should be make a hole in the surgical site through the drilling procedure.
- ③ To fixate plate and screw, the surgeon uses driver shafts. Insert the required driver shaft in the screwdriver body. Use it to pick up a screw then rotate it to fixate the screw in the surgical site.

### **3. Plate & Screw Removal**

In instances where Dual Top Anchor System removal is required, after normal surgical exposure the Dual Top Anchor System may be dislodged by using a driver shaft and other instruments. According to the surgeon's decision, if Dual Top Anchor System removal are required, they should be removed, inventoried and discarded. Do NOT reuse. The surgical site is now re-sutured.

### **4. Recommend usage site**

#### **① Maxilla**

- Inferior ridge of piriform aperture
- Posterior interdental septa
- Incisive interdental septa
- Inferior border of zygomatic process
- Posterior palatine alveolar bone
- Maxillary tuberosity
- Anterior rugae region
- Median palatine suture

#### **② Mandible**

- Alveolar bone
- Posterior Interdental septa
- Incisive interdental septa
- Retromolar area
- Body
- Parasymphysis

### **5. Orthodontic Appliance**

#### **① Dual Top Screw**

Type	Orthodontic Appliance	Driver shaft
JA series	Rubber band	MJ, MD Series
	Arch wire	
JB series	Rubber band	JB series
	NiTi Spring	
JD series	Arch wire	MJ series
G1 series	Arch wire	GD Series

G2 series	Rubber band	MJ Series
	Arch wire	
JF series	Rubber band	MJ, MD Series
	Arch wire	
JS series	Rubber band	MJ, MD Series
	Arch wire	

※ Orthodontic appliance is available when they are compatible with the product regardless of brand.

(Size : Rubber band – 1/8"~3/8", Arch wire – 0.012~0.018", 0.018X0.025"~0.021X0.025", Ni-Ti Spring : NCS0105, NCS0106, NCS0505) Only designated driver shaft is available for screw.

## ② Anchor Plate

Type	Orthodontic Appliance	Screw
AP-TL series	Rubber band	20-AT Series
AP-TR series		
AP-T series		
AP-L series		
AP-S series		

※ Orthodontic appliance is available when they are compatible with the product regardless of brand. (Size : Rubber band – 1/8"~3/8", Arch wire – 0.018X0.025"~0.021X0.025") Only designated screw is available for Plate.

## [INDICATIONS]

Dual Top Anchor System is used to provide a fixed anchorage point for attachment of orthodontic appliances to facilitate the orthodontic movement of teeth. It is used temporarily and is removed after orthodontic treatment has been completed. For use on adults over the age of 12.

## **[MATERIAL]**

Dual Top Screw : Titanium Alloy (ASTM F 136)

Anchor Plate : Titanium (ASTM F 67)

## **[WARNING]**

1. U.S. Federal Law restricts this device to sale by or on the order of a licensed surgeon or other qualified medical professional.
2. Selection of plate and screw size must be carefully considered by the operating surgeon and should take into consideration the quality of bone, bone type, functional loads exerted on bone(s), and post-operative patient compliance.
3. For safe and effective use of [Dual Top Anchor System], the surgeon should have specific training, experience, and thorough familiarity with the use of rigid fixation and techniques.
4. Multiple bending may weaken the plates and could result in implant fracture and failure.
5. The surgeon must exercise reasonable judgment when deciding which plate and screw type to use for specific indications.
6. Dual Top Anchor System are not intended to endure excessive abnormal functional stresses.
7. Use of an undersized plate or screw in areas of high functional stress may lead to implant fracture and failure.
8. Single use devices cannot be reused. Reuse of single use devices may lead to infection.

## **[CAUTION]**

1. Fixing the product incompletely may cause damage or breakage of the product, and the patient may feel a foreign body sensation.
2. Always follow appropriate safety precautions.
3. Select the appropriately sized plate and screw for the patient.
4. Responsibility for proper selection of patients, adequate training, experience in the choice and placement of plate & screw and the decision to leave or remove plate and screw postoperatively, rests with the surgeon.
5. Once applied, never reuse this device.
6. When engaging the screw, axial pressure of the driver shaft into the screw head must be adequately applied to ensure that the blade is fully inserted into the screw head. This results in proper axial alignment and full contact between driver shaft and screw, minimizing the risk of round-out. Otherwise there will be an increased risk of mechanical damage to the implant or the driver shaft.
7. Pilot drilling is required, if the bone density of the patient is too high or too thick.

8. Careful handling and storage of the product is required. Store in a dry place at room temperature.

#### [PRECAUTIONS]

1. The surgeon should inform the patient of risks associated with surgery.
2. Inspect each device to ensure they are not bent or damaged.
3. Before use, check that the specifications of the product matches those on the label.

#### [POSSIBLE ADVERSE EFFECTS]

1. In many cases, adverse results may be clinically related rather than implant related.
2. Osteoporosis, inhibited revascularization and poor bone formation can cause loosening, bending, cracking and fracturing of this device or premature loss of fixation with the bone, leading to nonunion.
3. Inadequate arrangements may cause procedure failure.
4. The product should not be used on patients sensitive to metals such as titanium, as it may cause boils and inflammatory reactions.
5. After inserting the product, there may be discomfort, pain, swelling, bleeding or root resorption on the patient.

#### [CONTRAINDICATION]

1. Not for use in cases of active or suspected infection or in patients previously sensitized to titanium.
2. Not for use in patients exhibiting disorders which would cause the patient to ignore the limitations of rigid fixation plate and screw implants.

#### [CLEANING]

Dual Top Anchor System's accessories (surgical instruments) can be reused after cleaning and sterilizing. They must be cleaned before reuse.

Following instructions are recommended for cleaning;

##### 1) Pre-cleaning

- ① Disassemble the device where possible.
- ② Remove gross soil using paper wipes and solution of cleaning agent.
- ③ Immerse the device in solution of cleaning agent for the time recommended by the detergent's manufacturer.

- ④ Using suitable brushes (never metal brushes or steel wool) cleaning the device thoroughly.
- ⑤ Rinse in running water until all traces of cleaning solution are removed.
- ⑥ Visually inspect for any remaining soil and repeat the steps above if necessary.

## 2) Cleaning

- ① Immerse the device completely in solution of cleaning agent and activate the ultrasonic bath at the concentration and temperature specified in the detergent manufacturer's instructions.
- ② Follow the same as the clause ④ ~ ⑥ of instructions for pre-cleaning.

## [STERILIZATION]

### • For Non-Sterile

Dual Top Anchor System must be sterilized before use. Following instructions are recommended;

- ① Place the devices in the appropriate block (kit or tray) using forceps and/or powder-free gloves to avoid contamination and any other negative effect on the surface of device.
- ② Wrap the block with a surgical drape.
- ③ Sterilize in the autoclave validated and maintained in accordance with ISO 17665 and ANSI AAMI ST79. Following parameters are validated in accordance with ISO 17665-1 and recommended for sterilization;

Cycle	Temperature	Exposure Time	Load Characteristics	Drying Time
Gravity	132°C	15 min.	Wrapped*	30 min.
Pre-vacuum	132°C	4 min.	Wrapped*	30 min.

\* In the case of load characteristics, we recommend usage of an FDA cleared wrap to ensure that the device is actually sterile prior to implantation.

- ④ Caution: Products intended for single use only must not be re-used or reprocessed. Devices labeled for single-use only should never be re-used or reprocessed. Re-use of these devices may cause harmful effects to patients. Re-use may compromise significant degradation in device performance, crossinfection, and contamination.

JEILMEDICAL Corporation does not take any responsibilities for any problems caused by re-use or reprocessing (e.g. cleaning and resterilization).

- **For Sterile**

Sterile Dual Top Anchor System is supplied by gamma sterilization.

Expiration date : Refer to label

#### [PACKAGING]

Dual Top Anchor System are packaged individually or set. This device manufactured and sold by JEILMEDICAL Corporation.

#### [SYMBOL DESCRIPTIONS]

<b>Symbol</b>	<b>Description</b>	<b>Symbol</b>	<b>Description</b>
	Catalogue number		Batch code
	Do not reuse		Do not resterilize
	Date of manufacture		Manufacturer
	Non-Sterile		Use-by Date
	Sterilized using irradiation		Caution
	Consult instruction for use		Do not use if package is damaged
	Authorized representative in the European Community		<b>EMERGO EUROPE</b> Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands

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Toto zariadenie vyrába a predáva spoločnosť JEILMEDICAL Corporation.