## INDICATIONS

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## INDICATIONS

These brackets are made from 99.9% pure alumina and are resistant to staining and discoloration. Their small size, smooth surfaces and rounded corners are designed for durability as well as patient's comfort. The use of the medical device is restricted to an orthodontic professional. Delivered in non-sterile state, the medical device is intended for single use.

## IDENTIFICATION OF THE BRACKET

These brackets are marked with colored dots for easy identification and placement. This color marking can be removed with alcohol. The color identifies the bracket/tooth combination. All colored dots appear on the distal gingival corner of the wing.

## STRAIGHT WIRE TECHNIQUE

Maxillary	Tooth	Torque	Angulation	Color
Incisor	1R	12/17	5	Orange
Incisor	1L	12/17	5	Orange
Incisor	2R	8	9	Green
Incisor	2L	8	9	Green
Canine	3R	0	11	Red
Canine	3L	0	11	Red
Premolar	4R & 5R	-7	0	Blue
Premolar	4R & 5R	-7	0	Blue

Mandibular	Tooth	Torque	Angulation	Color
Previous	1R & 2R	0	0	Violet
Previous	1L & 2L	0	0	Violet
Canine	3R	-11	5	Yellow
Canine	3L	-11	5	Yellow
Premolar	4R & 4L	-17	0	Brown
Premolar	5R & 5L	-22	0	Black

## CONTRAINDICATIONS

- Brackets should not be used on patients with poor oral hygiene or on patients who are unwilling to follow the treatment plan. Brackets should not be applied to compromised tooth surfaces or to teeth with existing root resorption.
- Occlusal bonding stresses should be avoided due to the hardness of the substrate and the risk of dental damage. Patients should avoid difficult and/or heavy oral insertions throughout the treatment.
- In the case of MRI tests (Magnetic Resonance Imaging), these medical devices must be removed depending on the area examined as they may produce artefacts.

## PRECAUTIONS FOR USE

- **Positioning of the bracket** : They have been designed for easy placement. Position the brackets in the middle third of the tooth with the slots parallel to the occlusal plane. Smaller models and rounded edges avoid problems during occlusion. The incisal edge of the bracket, parallel to the arch wire slot, can be used as an incisal/occlusal reference.
- **Bonding of the bracket** : Successful bonding can be achieved using a ceramic mechanical bonding system currently available on the market. These brackets have a specially designed rounded base to achieve maximum bonding surface between the bracket and the adhesive. The base design results in excellent bonding due in part to the increased surface area of the bracket to the adhesive. **DO NOT** touch the bases without gloves: particles from the skin may adhere and affect the strength of the adhesive. If a bracket should break but is still functional, it can be sanded with a diamond bur.
- **Removal of the bracket** : It does not require any special tools or techniques for debonding. Safe debonding can be achieved by using one of the commercially available ceramic bracket debonding instruments.
- **Ligature** : The wings have been designed with sufficient length to hold the bands securely in place. All ceramic brackets have been designed in a dual configuration, allowing for full or half bracket ligation. All brackets are designed for the strength required and are rounded for comfort.
- **Storage** : No special storage conditions.
- **Elimination** : Disposal of all orthodontic appliances must comply with current environmental regulations for waste.

## INDICATIONS

These brackets are made from medical grade polyurethane. Its small size, smooth surfaces and rounded corners are designed to be durable as well as comfortable for the patient. The use of this medical device is limited to an orthodontic professional. Delivered non-sterile, this device is intended for single use.

## IDENTIFICATION OF THE BRACKET

These brackets are marked with colored dots for easy identification and placement. This color marking can be removed with alcohol. The color identifies the bracket/tooth combination. All colored dots appear on the distal gingival corner of the wing.

## STRAIGHT WIRE TECHNIQUE

Maxillary	Tooth	Torque	Angulation	Color
Incisor	1R	12/17	5	Orange
Incisor	1L	12/17	5	Orange
Incisor	2R	8	9	Green
Incisor	2L	8	9	Green
Canine	3R	0	11	Red
Canine	3L	0	11	Red
Premolar	4R & 5R	-7	0	Blue
Premolar	4R & 5R	-7	0	Blue

Mandibular	Tooth	Torque	Angulation	Color
Previous	1R & 2R	0	0	Violet
Previous	1L & 2L	0	0	Violet
Canine	3R	-11	5	Yellow
Canine	3L	-11	5	Yellow
Premolar	4R & 4L	-17	0	Brown
Premolar	5R & 5L	-22	0	Black

## CONTRAINDICATIONS

- Brackets should not be used on patients with poor oral hygiene or on patients who are unwilling to follow the treatment plan. Brackets should not be applied to compromised tooth surfaces or to teeth with existing root resorption.
- Occlusal bonding stresses should be avoided due to the hardness of the substrate and the risk of dental damage. Patients should avoid difficult and/or heavy oral insertions throughout the treatment.
- In the case of MRI tests (Magnetic Resonance Imaging), these medical devices must be removed depending on the area examined as they may produce artefacts.

## PRECAUTIONS FOR USE

- **Positioning of the bracket** : They have been designed for easy placement. Position the brackets in the middle third of the tooth with the slots parallel to the occlusal plane. Smaller models and rounded edges avoid problems during occlusion. The incisal edge of the bracket, parallel to the archwire slot, can be used as an incisal/occlusal reference.
- **Bonding of the bracket** : A successful bond can be achieved using a mechanical bonding system currently available on the market. Follow the instructions enclosed with the bonding product you have chosen, observing the following general instructions:  
**DO NOT** touch the bases without gloves, particles from the skin may adhere and affect the strength of the adhesive.  
**DO NOT** clean brackets with water. They can be cleaned with alcohol if dried very carefully with air before gluing.  
**DO NOT** use a plastic activator with these brackets, it will reduce the strength of the adhesive.
- **Removal of the bracket** : They do not require special tools or techniques for debonding. Safe peeling can be achieved using one of the many bracket peeling instruments available on the market.
- **Ligature** : The wings have been designed with sufficient length to hold the brackets securely in place. All composite brackets have been designed for dual configuration, allowing for full or half bracket ligation. All brackets are designed for the strength required and are rounded for comfort.
- **Storage** : No special storage conditions.
- **Elimination** : Disposal of all orthodontic appliances must comply with current environmental regulations for waste.

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