Gebrauchsanweisung  
Instructions d'utilisation  
Instructions for use

# Tiger Gel clear

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**Dear doctor,** we are honoured that you have chosen this innovative new product from Tiger Dental GmbH.

To ensure the safe and effective use of this product, you must carefully read and adhere to this user manual.

The user manual covers all necessary safety instructions but cannot describe all possible circumstances. Therefore, you will also find video clips on our website [www.tigerdental.com](http://www.tigerdental.com) to enable an audio-visual description of the purpose of use, mentioned in this manual. Additionally, our customer service is available for any additional questions.

The continual development of our products indicates our company's commitment to innovation. Accordingly, after several clinical uses, we recommend you carefully review the accompanying user manual (this can be found online under the "Downloads" link at [tigerdental.com](http://tigerdental.com)).

#### 1. Distributed by

Tiger Dental GmbH, Allgäustrasse 3, 6912 Hörbranz, Austria

#### 2. Manufacturer

P.L Superior Dental Materials GmbH, 22607 Hamburg, Germany

#### 3. Product description

1. Tiger Gel is a light cure resin gel for use in dental practice and dental laboratories. It is used in orthodontics to coat all kind of metal parts in the mouth to make them smoother and enhance the experience and comfort for the patient.

2. It is easy and precise to apply and therefore sparingly and not time consuming. And the material is immediately positional stable.

3. Tiger Gel cures in the wavelength range between 400 and 500 nm and can be cured with commercially available light curing units.

#### 4. Indications

- Coating of metal parts (e.g. wires, bands, brackets) in the mouth
- Creation of bite blocks

#### 5. Application

Before applying the Tiger Gel to the corresponding metal parts or plastic surfaces, these must be cleaned thoroughly.

The Tiger Gel is not self-adhesive. This enables the user to easily remove it. If additional support is required, the surfaces to be treated must also be roughened.

Depending on the power of the light curing lamp, the Tiger Gel takes 20-180 seconds to harden.

Remove the syringe cap and lock one of the provided NeedleTips.  
NeedleTips are for single use only! After use replace the syringe cap.

Tiger Gel is applied directly from the syringe. To ensure complete curing seek for a constant layering (ca. 1-2 mm) while applying the material. Cure each layer separately.

Curing time depends on the used curing unit:

Kulzer Unix SX: 180 s

Kulzer Translux: 40 s

#### 6. Storage instructions

Do not store above 25°C (77°F). Protect from direct sunlight. Do not use after expiry date.

Close syringe immediately after use to avoid exposure to light.

#### 7. Safety Instructions

- Avoid contact of uncured material with eyes and skin.
- If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- On sensitive persons, sensitization to the material cannot be excluded. The material should not be used any more, if allergic reactions are observed.
- Do not use for patients with allergic reactions against acrylate.
- Do not breathe polishing dust; use suitable mouth protective device or aspiration!

#### 8. Warranty

P.L. Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture. P.L. Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and P.L. Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the P.L. Superior Dental Materials GmbH product.

#### 9. Limitation of Liability

Except where prohibited by law, P.L. Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

#### 10. Contents

1 pc. syringe Light Curing Resin 3ml

6 needle tips gauge 20

#### 11. Description of the used symbols of ÖNORM EN ISO 15223-1

 European Conformity

 Manufacturer

 Date of manufacture

 Lot code

 Item number

 Attention! Follow the instructions for use

 Expiry Date


 <sup>25 °C</sup> Storage Conditions - Upper Temperature limit

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 **0482** 

 Manufacturer:  
P.L Superior Dental Materials GmbH,  
22607 Hamburg, Germany

 Vertrieb:  
Tiger Dental GmbH, Belruptstrasse 59,  
6900 Bregenz, Austria  
[www.tigerdental.com](http://www.tigerdental.com)