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General remarks

1. Prior to first use

New instruments always require basic cleaning prior to first use / sterilization. All instruments must be cleaned, disinfected and sterilized before each use.

2. Validated methods for cleaning, disinfection and sterilization

Effective cleaning and disinfection constitute an absolute precondition for proper sterilization of the instruments. The user is responsible for the sterility of the instruments. Therefore, make sure only validated methods are used for cleaning, disinfection and sterilization. Sterilization equipment must also be serviced and checked at regular intervals. The validated parameters for the cleaning and sterilization cycles must also be checked regularly.

3. Legal requirements: Also comply with the legal requirements valid in your country and the hygiene regulations of the medical practice or hospital. Improper handling and care or use for purposes for which they are not intended may result in premature wearing out of instruments..

4. Personnel: Persons using these instruments should be familiar with use and handling of instruments, accessories and related devices.

5. Medical devices : Pliers and others orthodontics instruments that come into contact with mucosa or pathologically changed skin.

Inspection and functional check

It is very important to inspect each surgical/dental instrument for ruptures, cracks or malfunction prior to each use.

Above all, areas such as the cutting edges, tips, closures, locks and catches, as well as all moving parts, must be carefully checked.

Never use damaged instruments. Do not carry out repairs yourself. Service and repairs must be done by properly trained personnel only. If you have any questions in this regard, ask the manufacturer or your technical medical department.

Protection of personnel and third parties

Protective gloves that comply with the specifications of Directive 89/686/EEC must be worn when handling all used and contaminated instruments. Contaminated instruments must be disinfected at the earliest possible stage of reprocessing to maximize the safety of personnel handling contaminated instruments.

Steps in instrument reprocessing

Pre-treatment

Remove coarse soiling from instruments immediately after use (within max. 2 hours) before further processing in a tray/cassette system. Contaminated instruments must undergo pre-treatment within a maximum of 2 hours after use.

- Do not place in NaCl solutions (otherwise risk of pitting or stress corrosion cracking).
- Use only a released solution combining cleaning and disinfectant agents, (i.e. with DGHM or RKI approval or CE label) with no protein-fixing effect (be sure to following
- manufacturer recommendations for mixing).
- Avoid overfilling instrument sieves and wash trays.
- Always reprocess jointed instruments in open position.

Use only a soft brush, with long handle as needed, for manual removal of coarse soiling. Never use metal brushes or steel wool.

Remember that the disinfectants used in pre-treatment are for personal safety only and cannot replace the subsequent disinfection step.

Automatic cleaning

Equipment: Cleaning / disinfection device, cleaning agent

1. Place jointed instruments in the device with opened joints so that water can drain out of cannulas and blind holes.
2. Set cycle, wash for at least X minutes* and rinse for at least X minutes* (*see manufacturer's instructions).
3. When removing the instruments check cannulas, blind holes, etc. for visible soiling. Repeat cycle or clean manually if necessary.

Manual cleaning

Equipment: Cleaning agent, brush, running water

1. Rinse surface soiling off the instrument thoroughly.
2. Use a brush to apply cleaning solution to all surfaces. Be sure to clean jointed instruments in both the opened and closed positions.
3. Hold instrument under running water. The running water must flow through the cannulas and blind holes must be filled and emptied repeatedly.

Comment: Use a brush suitable for reaching the entire surface when cleaning blind holes.

Disinfection

Both neutral pH and alkaline cleaning agents can be used. Disinfectant solutions can be used according to the instructions on the label (see manufacturer information).

The automatic cleaning can be followed by thermal disinfection at 93° C / 200° F for 10 minutes (see manufacturer information for thermal disinfectant).

Use deionized water for the final rinse if possible. This will avoid spots, deposits and corrosion on the rinsed items.

Drying

For a drying phase reached as part of a cleaning / disinfection cycle the temperature should not exceed 93° C.

Maintenance

Apply a very small amount of high-quality silicone spray or medical white mineral oil to the joints and movable parts.

Sort out dull or damaged instruments. Check for cracks and damage. Check for proper functionality.

Control and functional test

Check jointed instruments for ease of movement (avoid excessive play). Check locking mechanisms (locking wheel) for functionality.

All instruments: Carry out visual inspection for damage and wear. Cutting edges should be free of notching and uniform. Check long, narrow instruments (in particular jointed instruments) for damage.

Packaging Individual:

Standardized packaging material can be used. The bag must be large enough for the instrument so the closure is not strained.

Sets:

Sort instruments into appropriated trays or place on all-purpose sterilization trays. Cutting edges must be protected. Package trays using a suitable method.

Sterilization

The sterilizer and sterilization method must comply with the valid and applicable standards and directives. EN 13060 Small Steam Sterilizers distinguishes 3 classes of autoclave:

Type B for packaged, massive, hollow and porous devices

Type N for unpackaged, massive instruments

Type S for devices as listed by small sterilizer manufacturers

- Gravitational method (with sufficient device drying).
- Steam sterilizer acc. to EN 13060 or EN 285.
- Validés conformément à la norme EN ISO/ANSI AAMI ISO 17665.
- Maximum sterilization temperature 134° C ; plus tolerance as per EN ISO/ANSI AAMI ISO 17665)
- Sterilization time (exposure time at sterilization temperature):20 min minimum at 121° C or 5 min at 134° C

IMPORTANT: Rapid sterilization leads to high instrument wear levels.

Additional information: When sterilizing several instruments in a sterilization cycle, the maximum sterilizer charge must not be exceeded (see manufacturer information).

Storage

Store instruments in dry rooms to avoid formation of condensation.

Exceptions

Chrome-plated instruments must not be placed in the thermal disinfectant or ultrasound bath.

Guarantee

All of our products are designed and manufactured to comply with the most stringent quality requirements. We assume no liability for products that have been altered compared to the original, used for purposes for which they were not intended or used improperly.

Products that are repaired or altered at any company other than the manufacturer are no longer authorized to bear the CE symbol.

IMPORTANT: All returned instruments must be disinfected and sterilized, with a document which prove sterilization and with sterilization date.

Recommended care products

- *ECI0150 Medical whiteoil instruments 50 ml.*
- *ECM1070 IP-Spray Instruments 400 ml.*

PROOF OF DECONTAMINATION

In compliance with the applicable statutory provisions and, above all, for the protection of our employees, we require a signed "Proof of decontamination" for each return. Please ensure that this "Proof of decontamination" is filled in and enclosed with each product return (complaint / repair / other reason for return) and package the product such that no risk of injury exists for our personnel in incoming goods inspection while unpacking.

Article no.: Article designation:

If possible provide information on

Delivery note N°: of [date] Batch N°:

Reason for return:

We hereby confirm that

- The medical device(s) enclosed with the letter has/have not come into contact with blood, tissue or other bodily substances/fluids and hygienic safety can be confirmed by signing.
- The medical device(s) enclosed with the letter has/have come into contact with blood, tissue or other bodily substances/fluids and have been disinfected, cleaned and sterilized according to the applicable hygiene requirements for medical devices and the manufacturer's specifications, which we confirm by signing.

Specifications on disinfection, cleaning and sterilization:

- Disinfection and cleaning have been performed manually
- Disinfection and cleaning have been performed mechanically
- Steam sterilization (3 min. at 134°C / 15 min. at 121 °C)
- Other method (please specify)

The medical device(s) enclosed with the letter could not be decontaminated!

Reason:

Stamp from the institution

Date

Signature (and full name)

IMPORTANT : Should we not receive this document or a comparable confirmation, we reserve the right to return the goods to the sender at their own expense.

E07Q05ProofDecontamination-EN-Rev00

SURVEY ON MAINTENANCE METHODES OF ORTHODONTIC INSTRUMENTS

The informations below are specific to procedures that were implementing in your office. They will also help us to give you an appropriate answer to the issue met. Thank you in advance for you cooperation,

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1 – <u>CLEANING</u> :		<input type="checkbox"/> YES	<input type="checkbox"/> NO
<input type="checkbox"/> Manual			
<input type="checkbox"/> Automatic	Material : -----		
	Time of cleaning: -----		
PRODUCT(S) : -----			
WATER :	<input type="checkbox"/> tap water	<input type="checkbox"/> demineralised water	<input type="checkbox"/> Other -----
RINSING :	<input type="checkbox"/> YES	<input type="checkbox"/> NO	
If yes, which kind of water ? : -----			
DRYING	YES	<input type="checkbox"/>	NO <input type="checkbox"/>
	Tissue	<input type="checkbox"/>	
	Compressed air	<input type="checkbox"/>	
	others	<input type="checkbox"/>	-----

2 – <u>DECONTAMINATION</u> :		<input type="checkbox"/> YES	<input type="checkbox"/> NO
<input type="checkbox"/> Manual			
<input type="checkbox"/> Automatic	Material : -----		
	Time of decontamination : -----		
PRODUCT(S) : -----			
WATER :	<input type="checkbox"/> tap water	<input type="checkbox"/> demineralised water	<input type="checkbox"/> Other-----
RINSING :	<input type="checkbox"/> YES	<input type="checkbox"/> NO	
If yes, which kind of water ? : -----			
DRYING	YES	<input type="checkbox"/>	NO <input type="checkbox"/>
	Tissue	<input type="checkbox"/>	
	Compressed air	<input type="checkbox"/>	
	others	<input type="checkbox"/>	-----

3 – <u>STERILIZATION</u> :		<input type="checkbox"/> YES	<input type="checkbox"/> NO
<input type="checkbox"/> Cold sterilization – product :	-----		
<input type="checkbox"/> Autoclave	Time of cycle : -----	/ T°maximum : -----	
	Protection : <input type="checkbox"/> YES – which way (bags, oiling...) :	-----	
	<input type="checkbox"/> NO		

4 – <u>LUBRIFICATION</u> :		<input type="checkbox"/> YES	<input type="checkbox"/> NO
When ? -----	/ Frequency ? -----		
Product(s) : -----			

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