9

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Then heat the BIOPLAST® 2,0 mm by setting the code or recommended heating time, lock the pressure chamber and open after cooling phase has elapsed. Cut plate at model bottom.

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Heat the occlusal surfaces under the infrared heater using the heating support with the pressure chamber closed, but not locked. Insulate occlusally the antagonist model with BIOPLAST® insulating agent, in order to avoid plaster residues. Then put models in the articulator for counter bite. Cool down with compressed air and water and remove.

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Cut the plate with scissors to desired extensions.

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Polish the edges with BIOPERM trimmer or DIMO®PRO wheel. Cut the mouthguard diagonally to remove sharp edges in the labial and palatal area. Finish the rims using the OSAMU polisher at a low RPM exerting low pressure.

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Innovation for safety plus.

The new Mouthguard type BIOPLAST® XTREME gives safety plus in sports thanks to the integrated reinforcement in the delicate anterior region. Hard and flat impacts can be absorbed in contact sports, using this innovative new plate. The innovation: For the first time, different degrees of hardness are integrated in one pressure moulding plate. The size of the anterior protection area can be varied individually thanks to the V-shaped hard insert.

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Three Mouthguard types – meeting all requirements.

BIOPLAST® MASTER:
• proven mouthguard of two layers for maximum safety
• for all leisure time activities

The product innovations BIOPLAST® XTREME and BIOPLAST® XTREME PRO give safety plus thanks to specific extras.

BIOPLAST® XTREME:
• with integrated reinforcement segment in the anterior region
• for contact sports; absorbs hard and flat impacts

BIOPLAST® XTREME PRO:
• with integrated reinforcement segment plus hard-elastic layer in the anterior region
• for professional contact sports where hard and flat impacts often occur, absorbing the impact power over a large surface.

All mouthguard types as described above consist at least of 2 layers, guaranteeing maximum safety and comfort.

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Apply one thin layer of BIOPLAST® insulating agent only to those parts of the model which will be covered by the mouthguard. Trim model flat and place on model platform with the clear V-shaped part of the BIOPLAST® XTREME PRO plate covering the desired extension of the anterior section (e.g. from 3-3).
 Make sure the side of the label with the imprint is facing the infrared heater when fixing the plate. Heat the foil by setting the code or recommended heating time according to the instructions. Close pressure chamber and open after cooling phase has elapsed.

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Cut plate with scissors or scalpel on the model bottom and insulate the anterior section with a silicone spray.

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For the fabrication of a protection insert: embed the model vertically in the pellets and heat a plate DURAN® 0,75 mm by setting the code or recommended heating time. Lock pressure chamber and open after cooling phase has elapsed.

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Cut the BIOPLAST® XTREME to the desired extensions. Make sure the plate is being cut straight edge, as cutting it diagonally might impair the optimum lamination with the 2nd BIOPLAST® plate.
 Cut the DURAN® plate, too, with scissors to the desired extensions.

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Make sure the extension of the DURAN® plate towards the alveolar ridge is some millimeters shorter than the extension of the BIOPLAST® XTREME PRO underneath, in order to avoid problems in laminating the plates.

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Pierce holes into the BIOPLAST® XTREME PRO plate - occlusally and interdentially - by means of a bur, in order to avoid air bubbles under the protection insert and between the BIOPLAST® plates.

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Clean the BIOPLAST® XTREME PRO plate with moist tissue paper to remove the silicone insulation.

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Place the model covered with the protection insert for the anterior section made of DURAN® with the front in the direction of the pressure chamber.

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