**FOR PROFESSIONAL USE ONLY**

Instructions for Use

**INDICATIONS FOR USE:**

FLI® Ceramic Adhesive is designed for use as an adhesive for retaining orthodontic brackets onto tooth enamel.

**DESCRIPTION:**

FLI® Ceramic Adhesive is a light cured orthodontic bracket adhesive system. It includes a thick viscosity composite adhesive, a low viscosity resin sealant and an etchant gel. The light curing mechanism assures a quick setting to a strong adhesive system.

**FORMULATION:**

FLI® Ceramic Adhesive is a composite based on a BisGMA resin system with a small particle glass filler and a photosensitive catalyst system. FLI® Ceramic Adhesive Sealant is an unfilled resin system based on BisGMA with a photosensitive catalyst system.

The etchant is a 37% phosphoric acid gel.

**INSTRUCTIONS:**

**PREPARATION**

- Prophyl the teeth with a non-fluoride, oil-free prophyl paste. Wash with an oil-free air/water spray. Do not allow the patient to rinse.
- Isolate teeth with rubber dam or cotton rolls.
- Blow dry completely with an oil/moisture-free air source.

**ETCHING**

- Dispense 1 or 2 drops of the etchant into a mixing well. One may also apply the etchant directly to the tooth with a syringe and a syringe tip. Recap the etchant.
- Apply the etchant to the enamel with a brush or sponge pellet. Dab the etchant onto the tooth, but do not rub against the enamel. Allow etchant to remain for 30 seconds.
- With a few drops of water from the air/water syringe and with high volume suction near the tooth, wash away the gel. This avoids splashing the etchant. After the etchant is gone continue to wash the area of enamel completely for 20 seconds with an air/water spray and adequate suction. Do not allow saliva to contaminate etched enamel surfaces. Contaminated surfaces need to be re-etched for 10 seconds.
- Air dry thoroughly using oil/moisture-free air. After drying the etched surfaces of teeth should have a dull, white, frosted appearance. If the enamel is not chalky appearing, repeat the process allowing etchant to remain on the teeth for additional 30 second increments of time and wash and dry as described in this etching procedure. This repetition may be needed especially on teeth high in fluoride content.

**SEALING**

- Dispense 2 drops of Sealant liquid into a clean mixing well. Recap the Sealant bottle.
- Apply the Sealant directly onto each etched, dry surface. Remove excess Sealant Resin with a brush or a gentle stream of oil-free, moisture-free air. Avoid contact with gingiva or other soft tissue.
- Curing the Sealant is not necessary at this time but can be done to limit the mobility of the bracket. Hold the end of the light guide as close to the tooth as possible without touching the resin. Expose the Sealant for 10 seconds on each tooth. (WARNING: Use a blue curing light, bandwidth of 470-480 nanometers with a minimum of 400 mW/cm<sup>2</sup> intensity.)

INST 604-0

**ADHERING**

- Apply a small increment of Adhesive Paste directly onto the bracket back, or use a small, flat ended hand instrument to place the composite on the bracket. Make sure the paste is worked into the bracket backing. Recap the syringe or syringe tip after use.
- Seat the bracket onto the tooth with light pressure and position the bracket. Press firmly to place but do not express all of the paste from under the bracket.
- Remove any gross excess of paste that expresses from around the bracket.
- Light cure the Adhesive from two directions (incisal and gingival or mesial and distal) for 20 seconds each direction for each bracket. The light guide tip should be within 5 mm of the tooth surface. (WARNING: Use a blue curing light, bandwidth of 470-480 nanometers with a minimum of 400 mW/cm<sup>2</sup> intensity.)
- Arch wires can be placed immediately after bonding the last bracket.

**BRACKET REMOVAL**

- At completion of the orthodontic treatment remove the bracket by placing the cutter edges at the adhesive/tooth interface. Remaining adhesive can be removed with a rubber wheel and/or scaler.

**WARNINGS:**


- Dentists and assistants should wear gloves and protective eye wear. Patients also should wear eye protection.
- Etchant gels are caustic and can cause tissue damage in contact with skin or eyes. If accidental contact occurs with skin, flush immediately with copious amounts of water. For eye contact, flush immediately and for a prolonged period with copious amounts of water and consult a physician.
- Uncured monomers in this resin material may cause an allergic reaction or skin irritation in some individuals. Do not use in individuals with allergies to these materials.
- Any non-sterilizable item that is handled in the delivery of the dental service, should be disinfected by standard dental office hygiene procedures.
- Do not allow saliva or water contamination of etched tooth surface or unset material.
- If a fluoride treatment is planned at the same appointment, apply the fluoride after bonding the bracket, and the adhesive is completely set.

**STORAGE CONDITIONS:**


The system is designed to be stored between 10°-24° C (50° -75° F). Replace caps immediately after dispensing. If stored in a refrigerator, bring to room temperature prior to use. Do not store materials in the proximity of eugenol-containing products. Do not store the materials in intense light, or under wet conditions.

*Manufactured for*  
Rocky Mountain Orthodontics™  
650 West Colfax Avenue  
Denver, CO 80204 USA  
303.592.8200/  
Toll Free 800.525.6375

Made in the USA  
Fabricado en los EE.UU.  
Fabriqué aux Etats-Unis  
Hergestellt in den USA

European Representative  
RMO Europe  
Rue Geiler de Kaysersberg  
67400 Illkirch, France

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