
## Set de finition pour la technique de thermoformage Finishing-Set for pressure moulding materials

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Art.-Nr. 3247

polissage des plaques sandwich

### Finishing-Set for pressure moulding materials

Art.-No. 3378

#### Cutting bur (white)

Art.-No. 3214

Rough cutting of hard foils

#### HM-carbide cutter (yellow)

Art.-No. 3369

Rough cutting of hard foils and hard / soft foils

#### HM-Tricutter (green)

Art.-No. 3370

Finishing of hard foils and splints

#### Finishing bur (red)

Art.-No. 3377

Cutting of hard foils and splints with integrated wire elements.

#### Bioperm-Trimmer (blue)

Art.-No. 3226

Rough trimming of soft foils and hard / soft foils

#### Polyfix (black)

Art.-No. 3371

Finishing of soft foils and hard / soft foils

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IMPRELON® »S« in the a. m. thicknesses is a high quality polycarbonate which is a break- and abrasion-resistant clear thermoplastic. The material is best suited for all splint therapies.

**Pre-drying:** Due to the properties of polycarbonate and the high pressurizing temperature of approximately 180 °C, **pre-drying** of IMPRELON® »S« is absolutely **necessary** to **avoid air bubbles** when plastifying this material. IMPRELON® »S« can be stored in the BIOSTAR-drying oven or, depending on the thickness of the material, can be **pre-dried** for the recommended conditioning time in the BIOSTAR drying oven or a similar oven at a temperature of **80 °C – 100 °C**. It may be necessary to prolong first pre-drying time of 2.0 and 3.0 mm sheets depending on humidity, storing and transport conditions. Sheets prepared in this way can be used for approximately 7 days without any air bubbles. Repeated pre-drying is possible without any loss of quality.

**Bite splints:** The etching liquid PERLIBOND can be applied to IMPRELON® »S« if resin should be added to bite splint. For this use 1 – 2 times PERLIBOND for the occlusal areas and model the desired shape by applying autopolymerisate (e. g. DURASPLINT®, STEADY-RESIN) on the moistened bonding agent.

Polymerization and polishing is done as usual.

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IMPRELON® »S« – as well as all our pressure moulding foils – has been tested at the University of Heidelberg with regard to toxic properties, gingiva irritations and mutagenity. The results clearly state that these foils are biologically well-accepted.

## Programme de livraison / Ordering

Art.-No.	IMPRELON® »S«	Temps de présechage recommandé Recommended Pre-drying time
<b>3357</b>	0,50 x 125 mm O	1.5 h
<b>3171</b>	0,75 x 125 mm O	2.0 h
<b>3172</b>	1,00 x 125 mm O	4.0 h
<b>3389</b>	1,50 x 125 mm O	5.0 h
<b>3276</b>	2,00 x 125 mm O	7.0 h
<b>3275</b>	3,00 x 125 mm O	12.0 h
<b>4239</b>	<b>PERLIBOND</b>	
<b>3129</b>	<b>BIOSTAR®- / MINISTAR®-four de séchage / Drying Oven</b>	