


# ALPINE SL™ CLEAR

Fully interactive aesthetic self-ligating bracket, with Nickel-Titanium spring clip



## FEATURES

ADVANCED CERAMIC MATERIAL

NITI CLIP WITH RHODIUM-COATING

EASY OPEN/CLOSE FUNCTION

FULL SLOT-WIDTH  
ENHANCED CLIP DESIGN

SMOOTH, LOW-PROFILE SURFACE  
TIE-WING UNDERCUTS

INK MARK

CONTOURED BASE

TRIPLE-X MECHANICAL BASE

## BENEFITS

- Strong, fracture-resistant material with an aesthetic appearance
- Shape alloy ensures reliable clip operation and control throughout entire treatment; Rhodium-coating for aesthetically, inconspicuous appearance
- No special instrument required; Closes gingivally with slight pressure using figure on the incisal curve or using an instrument
- Clip fully covers slot for optimal rotational control
- Three-finger design secures archwire in place while preventing mesial-distal movement; Clip fully extends through the bracket for overall structural integrity
- Offers exceptional patient comfort
- Allows optional ligation with chains/ties
- Easy bracket identification
- Anatomical fit for greater control
- Triple-X base, multi-contoured design with undercuts to ensure precise fit and secure mechanical-bond

# ALPINE SL™ CLEAR

Fully interactive aesthetic self-ligating bracket, with Nickel-Titanium spring clip



## ALPINE SL™ CLEAR • ROTH PRESCRIPTION

MX	Torque	Angulation	Offset	Hook	R/L	.022
Central	12°	5°	0°		UR	Q12400
	12°	5°	0°		UL	Q12401
Lateral	8°	9°	0°		UR	Q12402
	8°	9°	0°		UL	Q12403
Cuspid	0°	11°	4° Mesial	D	UR	Q12406
	0°	11°	4° Mesial	D	UL	Q12407
Bicuspid	-7°	0°	0°	D	UR	Q12410
	-7°	0°	0°	D	UL	Q12411
MN	Torque	Angulation	Offset	Hook	R/L	.022
Anteriors	0°	0°	0°		LR/LL	Q12420
Cuspid	-11°	5°	2° Mesial	D	LR	Q12426
	-11°	5°	2° Mesial	D	LL	Q12427
<b>Description</b>						<b>.022</b>
Alpine SL Clear Roth 022 - Kit 5-5 With Hook 3/4/5						<b>ECKQ12402</b>

## ALPINE SL™ CLEAR • MBT PRESCRIPTION

MX	Torque	Angulation	Offset	Hook	R/L	.022
Central	17°	4°	0°		UR	Q12100
	17°	4°	0°		UL	Q12101
Lateral	10°	8°	0°		UR	Q12102
	10°	8°	0°		UL	Q12103
Cuspid	-7°	8°	0°	D	UR	Q12106
	-7°	8°	0°	D	UL	Q12107
Bicuspid	-7°	0°	0°	D	UR	Q12110
	-7°	0°	0°	D	UL	Q12111
MN	Torque	Angulation	Offset	Hook	R/L	.022
Anteriors	-6°	0°	0°		LR/LL	Q12120
Cuspid	-6°	3°	0°	D	LR	Q12126
	-6°	3°	0°	D	LL	Q12127
<b>Description</b>						<b>.022</b>
Alpine SL Clear MBT 022 - Kit 5-5 With Hook 3/4/5						<b>ECKQ12102</b>