

NOTICE D'INSTRUCTION

Ce dispositif médical est à usage unique.

Le risque d'une contamination croisée existe si une décontamination /stérilisation conforme à la réglementation n'a pas été réalisée.

Précaution.

Ce dispositif médical ne doit pas être utilisé sur un patient connu comme étant allergique au nickel.

Nettoyage.

La bague ainsi que les attachements est livrée propre et non stérile.

Cependant avant insertion en bouche nous recommandons son nettoyage par ultrasons.

INSTRUCTION FOR USE***This Medical device is for single use***

A cross contamination can occur if reused and if a decontamination/sterilization in accordance with regulation has not been realized.

Caution

This medical device contains nickel and should not be used for individuals with known nickel sensitivity.

Cleaning

Band and its attachments are delivered clean and non-sterile. However before inserting it into the mouth, we recommend its ultrasonic cleaning.

Description des symboles / Symbols description



Référence catalogue
Catalogue number



Code lot
Batch code



Fabricant
Manufacturer



Ne pas réutiliser
Do not reuse



Contient Nickel
Contains Nickel



Consultez les instructions d'utilisation
Consult instruction for use



Non stérile
Non-sterile